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“English Version”

Introduction:

The cleaning efficacy of two wash detergents were performed using the protocol for wash detergents V020712 already agreed with Ruhof.

The study requested was to perform a quantitative surgical instrument study in an SSD.

Rat brain tissue was used in the study as it has been demonstrated in our laboratory that proteins in brain tissue strongly adheres to stainless steel and provides a realistic test of the cleaning efficacy of wash detergents. Single use stainless steel instruments were used for the study. Analysis used the novel fluorescent protein detection technology developed at Barts where the method is optimised and validated.

Materials

Rat brains: Rat brains were collected from male Wistar rats (Charles River Labs). OPA/NAC reagent: Patented reagent prepared in the lab ProReveal System: Patented protein detection System. . Surgical Instruments: Purchased from Surgical Holdings Ltd (U.K). Washer: BHT Hygienetechnik Wash detergents:

- 1) Endozime XP (Ruhof)
- 2) Major Alkaline brand

Wash cycle

Programme: See appendix A for full details.

Cleaning time in the presence of Endozime XP = 10 minutes 50 seconds. Temperature: 45°C – 55°C Dilution of detergent

- 1) Endozime XP : 2ml/L in R0 water
- 2) Alkaline: 4ml/L in R0 water 2 R130812 (Dr.N.K.Nayuni)

Method:

Quantitative surgical instrument study (In an SSD wash using BHT Hygienetechnik washer disinfector) 40 new single use surgical instruments purchased for the study.

The instruments were initially washed in an SSD to remove any manufacturing debris and potential protein contamination due to handling. The washed instruments were packed into sterile bags in the controlled clean room environment of an SSD before transport to the laboratory. The instruments were contaminated with rat brain tissue by manipulating whole brains. The weight of adhered tissue was determined using a 4-place analytical balance. The brains were then homogenised and the total protein measured using a Biuret reagent (Bio-RAd Ltd). Following the manipulations the instruments were allowed to dry at room temperature for 48 h. The instruments were immediately transported to the SSD where they were cleaned in one of two hospital washer disinfectors (BHT Hygienetechnik) using standard cycles.

Each AWD had been optimised for one of the wash detergents by engineers using tests according to the manufacturers' guidelines. The washed and dried instruments were sealed in sterile bags in the clean SSD environment and returned to the laboratory at Barts for analysis. The residual protein on each

instrument was measured against a protein (BSA) standard using the OPA/NAC protein detection method.

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